



GENERATOR CONSOLES


Mini-Console



Membrane Rad Console






Membrane R&F Console



Dimensions:
mm (inches)

GENERATOR CONSOLES

StarQ Touch Screen Console (15.6"/400 mm)



(Additional specifications available on request)

Dimensions:
mm (inches)

CPI Console/Generator Compatibility

CONSOLE	Generator		
	CMP 200® SERIES	CMP 200®DR SERIES	INDICO IQ® R&F
Mini Controller	Not Available	Available	Available
Membrane Rad	Available	Available	Not Available
Membrane R&F	Not Available	Not Available	Available
StarQ Touch Screen (15.6"/400 mm)	Not Available	Available	Available

CPI Console Options

CONSOLE	Options			
	HAND SWITCH	WALL MOUNT	DESKTOP BASE	FLOOR PEDESTAL
Mini Controller	Included	Not Available	Included	Not Available
Membrane Rad	Available	Available	Not Available	Available
Membrane RF	Available	Available	Not Available	Available
StarQ Touch Screen (15.6"/400mm)	Included**	Available	Included	Not Available

** Includes mini controller & hand switch



Registered Trademarks are U.S. Registered Trademarks of Communications & Power Industries Canada Inc.

RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
 613-799-4171

