



CMP 200[®] ***SERIES***

High Frequency X-Ray Generator

*Leading Technologies.
Powered by*





Leading technologies have been powered by CPI performance and reliability reaching back to 1948. Through ongoing product development, support and customer service, CPI leads the field globally in innovation, providing advanced features, integration, flexibility and customization.

The CMP 200[®] High Frequency X-Ray Generator Series continues CPI's history of innovation by providing a high performance solution with enhanced features for the radiographic market that is unmatched in reliability.






BENEFITS

- Suitable for Film, CR and DR applications
- Advanced control interfaces for Digital Radiography
- 1 Φ and 3 Φ mains input power options
- Energy Assist capacitor storage module available

PERFORMANCE

- Wide output power range from 32 kW to 80 kW
- Mains input power ranging from 1 Φ /208-230 VAC to 3 Φ /400-480 VAC
- kV-AEC, kV-mAs, kV-mA-ms or kV-cm thickness techniques
- Exposure output from 125 kV to 150 kV, 10 mA to 1000 mA, 0.1 mAs to 1000 mAs
- Exposure time from 1 ms to 30 seconds
- Advanced options including AEC, DAP, Dual Energy, Serial Radiography, Stitching available
- Compatible with over 300 x-ray tube models using Low or Dual Speed Starter

Consoles Available



RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
613-799-4171

