



FILM/CR IMAGING SOLUTIONS

CH1 17" X 17" UPRIGHT IMAGING

Budget Friendly, No Compromise in Quality

TXR offers one of the most affordable upright systems that incorporates safety, reliability, value, quality and ease of use. At TXR we recognize that your time is valuable. The TXR CH1 system includes features that will increase patient thru-put and minimize the time needed for exams. Experience the TXR difference with standard features that other manufacturers option.

STANDARD FEATURES INCLUDE:

Floor Wall Mount Tube Stand

- Platform tube mount, 90° tube rotation
- Magnetic locks
- Angle indicator
- 72" Floor Track

Wall Stand

- Floor wall mounted
- Magnetic locks
- 17" x 17" Grid cabinet
- 103 line, 10:1 ratio grid
- Cassette tray

Generators

- 320, 400 or 500 mA High frequency



Collimator

- Certified manual collimator
- 25 Second timer
- Laser light
- Retractable tape measure
- Swivel mount

Options

- Super speed bucky
- Various grids
- Tilt mechanism
- Head clamp
- Side mounted grid cabinet or bucky
- Fixed SID, see back page




FILM/CR IMAGING SOLUTIONS

OPTIONAL FIXED SID

Budget Friendly, No Compromise in Quality

TXR offers one of the most affordable upright systems that incorporates safety, reliability, value, quality and ease of use. At TXR we recognize that your time is valuable. The TXR CH1 system includes features that will increase patient thru-put and minimize the time needed for exams. Experience the TXR difference with standard features that other manufacturers option.



STANDARD FEATURES INCLUDE:

Fixed Tube Stand

- Platform tube mount, 90° tube rotation
- Magnetic locks
- Angle indicator
- Mounts floor to wall fixed at 40" or 72" SID

Wall Stand

- Floor wall mounted
- Magnetic locks
- 17" x 17" Grid cabinet
- 103 line, 10:1 ratio grid
- Cassette tray

Generators

- 320, 400 or 500 mA High frequency

Collimator

- Certified manual collimator
- 25 Second timer
- Laser light
- Retractable tape measure
- Swivel mount

Options

- Super speed bucky
- Various grids
- Tilt mechanism
- Head clamp
- Side mounted grid cabinet or bucky



RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
613-799-4171