

# Weight Bearing Protection

*Our new Weight Bearing Protective Cover has been designed to accommodate Digital DR Panels and CR Cassettes.*

*Our covers are strong yet lightweight, offering a high degree of protection during weight bearing procedures.*

*Our latest cover is constructed of a high density polyethylene frame and a carbon fiber top, offering protection up to 750 lbs. Our easy lift handle makes removal and carrying a breeze.*


*The Weight Bearing Cover has been tested by an accredited independent laboratory.*



# Weight Bearing Protection



Polycarbonate



Carbon Fiber

- Available in Polycarbonate or Carbon Fiber
- Rated up to 750 lbs.
- Easy lift handle
- Non-skid rubber feet
- Multiple sizes available
- Printed target area

***This product is perfect for bariatric, equine, and podiatry procedures.***



Easy Lift Handle

Protect your **investment** use a **cover!**



**RC Imaging**  
your radiology resource

RAYTECH DIAGNOSTICS  
<https://www.raytechdiagnostics.com/>  
[sales@raytechdiagnostics.com](mailto:sales@raytechdiagnostics.com)  
613-799-4171

