



RAYMAX

MEDICAL IMAGING



Raymax introduces an affordable Podiatric Digital Radiography (DR) solution with superior quality. The Raymax Podiatric DR unit can be integrated with your existing X-Ray equipment and compatible with existing patient software.



Specifications

- ✓ 14x17 inch field of view in regular mode.
- ✓ 12 mega-pixel resolution (3024x2016 true pixels).
- ✓ .006 inch pixel dimension, 7 pixels per millimeter.
- ✓ 4 LPM
- ✓ .004 inch pixel dimension, 10 pixels per millimeter.
- ✓ Imager dimensions 15"x19"x22"

Computer Specifications

- ✓ Computer (minimum 1 Gig of RAM and minimum 250 Gig of hard disk storage).
- ✓ CD & DVD burner.
- ✓ 20 inch rotating LCD monitor with 1680x1050 resolution.
- ✓ Remote foot switch to capture X-Ray.
- ✓ UPS battery back-up.
- ✓ All required cables.
- ✓ Completely configured with DR

X-rays are saved in jpeg, bmp, or dicom format for easy importing into any office management or patient software on any computer.



X-rays may be transferred to any computer system via disk, network, or internet. X-rays may be transferred to a tablet pc with the DR software to adjust the x-rays with a touch of the tablet pen.



The high resolution images may be transferred to all types of media such as projectors or plasma televisions for patient viewing or presentations. Requires video out capability from a computer system.



DISTRIBUTED BY:

RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
 613-721-1982