

TigerView PACS-NCA

Acquire | Store | View | Share

TigerView PACS Software for NUCCA



TigerView allows you to uniquely integrate and manage digital x-rays in to your workflow leading to efficient patient education and patient care.

- *DICOM 3.0, HL7 interface, modality worklist and multiple viewing stations*
- *Unmatched user-friendly interface*
- *Multiple hardware manufacturer integration*
- *Multiple modality support*
- *Multiple exam room viewing*
- *NUCCA Analysis specific tools and measurements*
- *CD burning software with built-in CD viewer*
- *Integrates with leading Practice Management Systems*
- *Cloud viewing on PC's and iPad/iPhone's along with Cloud backup*
- *Customizable report generation*
- *Complete turn-key solution: set-up, staff training and operation of digital x-ray has never been easier*

Call us for an in-office demo today.

1-800-385-9593

TigerView PACS-NCA

TigerView PACS Software for NUCCA

SPECIFICATIONS

TigerView Functionality

Key DICOM 3.0 services, including: DICOM Receive Storage Commit and Modality Worklist
DICOM send
View DICOM header information
DICOM CD/DVD import and archival
DICOM CD/DVD burn functionality
Ability to edit patient/study information
Easy to filter, search and sort modality worklist
Worklist search
Lossless image compression
Multiple hardware manufacturer integration
Customizable Report generation
View images from any web-enabled PC w/TigerView Cloud
Touchscreen friendly PAC
Dual-monitor friendly
User customizable filters
Multiple image viewing

PACS Viewer - Modality Support

Direct Digital Radiography (DR)
Computed Radiography (CR)
Digitized x-ray film
Ultrasound (US)
Intra-oral cameras
Extra-oral cameras
Color digital camera images for morphing & wound care progression

EMR or Charting Bridging - Easy link from EMR

Proprietary bridges
DICOM modality worklists
HL7
TigerLink

Full set of editable NUCCA measurement tools

Perpendicular Baseline Measurement	Image Marking and Text Box Features
Baseline Measurement	Mark Spots
NUCCA Analysis	Center Point
Torgs Ratio	Stitching
Illum Analysis	Automatic Calibration Tool
Jacksons Angle	Cropping
L/R Cervical Normal Arc	Magnification and Zooming
L/R Thoracic Normal Arc	Window Leveling
L/R Lumbar Normal Arc	Colorized Density
Georges Line	Image Negation
Cobbs Angle	Image Sharpening
Vertical Offset	Spine Labeling
Horizontal Offset	AC - Atlas Circles
CR - Cranial Rotation	PL - Plane Line
CSL - Central Skull Line	N-T - Nasium Tracker
APL - Atlas Plane Line	V-T - Vertex Tracker
CC - Condyle Circles	

TigerView Cloud

Remote viewing
Multi-office and multi-provider sync
Centralized database
Disaster recovery
Off-site Cloud storage/backup
Automatic software updates

For complete product information, visit televeresystems.com