



- DICOM 3.0, HL7 interface, modality worklist and multiple viewing stations
- Unmatched user-friendly interface
- AVDF winning dental solution multiple years
- Supports full body and dental radiographs and photographs all in one platform
- Multiple hardware manufacturer integration
- Multiple exam room viewing
- Veterinary specific tools and measurements including vertebral heart score and hip dysplasia measurements
- Integrates with leading Veterinary Practice Management Systems
- Cloud viewing on PC's and iPad/iPhone's along with Cloud backup
- CD burning software with built-in CD viewer
- Customizable report generation
- Complete turn-key solution: set-up, staff training, and operation of digital x-ray has never been easier

TigerView PACS Software for Veterinarians

TigerView allows you to uniquely manage full body and all dental images seamlessly with one database



SPECIFICATIONS

TigerView Functionality

Key DICOM 3.0 services, including: DICOM Receive Storage Commit, and Modality Worklist

DICOM send

View DICOM header information

DICOM CD/DVD import and archival

DICOM CD/DVD burn functionality

Ability to edit patient/study information

Easy to filter, search and sort modality worklist

Worklist search

Lossless image compression

Multiple hardware manufacture integration

Full set of editable Measurement & Annotation tools:

- Vertebral Heart Score
- Cobb's Angle
- Hip Dysplasia
- Automatic Calibration Tool
- Cropping
- Magnification and Zooming
- Window Leveling
- Colorized Density
- Image Negation
- Image Sharpening
- Image Marking and Text Box Features

Customizable Report generation

View images from any web-enabled PC with TigerView Cloud

Touchscreen friendly PAC

Dual-monitor friendly

User customizable filters

Multiple image viewing

PACS Viewer – Modality Support

Direct Digital Radiography (DR)

Intra-oral cameras

Computed Radiography (CR)

Extra-oral cameras

Digitized x-ray film EMR

Color digital camera images for morphing and wound care progression

Ultrasound (US)

EMR or Charting Bridging – Easy link from EMR

Proprietary bridges

HL7

DICOM modality worklists

TigerLink

TigerView Cloud

Remote viewing

Disaster recovery

Multi-office and multi-provider sync

Off-site Cloud storage/backup

Centralized database

Automatic software updates

For complete product information, visit televersystems.com.

RAYTECHDIAGNOSTICS
<https://www.raytechdiagnostics.com>
sales@raytechdiagnostics.com
 613-799-4171