

# TigerView Pro PACS for Dentistry

**TigerView**  
SOFTWARE

Acquire | Store | View | Share



**TigerViewCloud**  
Upload | Store | View | Share



TigerView offers a complete digital workflow solution for image capture from almost all devices, image retrieval and viewing.

TigerView is a simple yet powerful imaging acquisition software to control and acquire images from all of the popular imaging devices that are used in the dental market.

The TigerView software is a complete image management system allowing you to acquire, store, view and share images and provides functionality that will significantly improve your work flow with better patient management and accurate diagnoses. Quickly and easily access all digital image files and share them with colleagues and patients. Image manipulation and specialty specific enhancement tools provide maximum information for minimum effort.

- DICOM 3.0, HL7 interface, modality worklist and multiple viewing stations
- Unmatched user-friendly interface
- One database for all imaging products
- Supports all sensors, digital pans, intro-oral cameras, film scanners all in one platform
- Multiple hardware manufacturer integration
- Multiple exam room viewing
- Dental specific tools and measurements
- Integrates with leading Dental Practice management systems
- Cloud viewing on PC's and iPad/iPhone's along with Cloud backup
- CD burning with built-in CD viewer
- Customizable report generation
- Complete turn-key solution: set-up, staff training, and operation of digital x-ray has never been easier

*Call us for an in-office demo today.*

**1-800-385-9593**

**TELEVERE**  
SYSTEMS



# TigerView Pro PACS for Dentistry

TigerView allows you to uniquely integrate and manage digital x-rays in to your workflow leading to efficient patient education and patient care.

## SPECIFICATIONS

### TigerView Functionality

- Key DICOM 3.0 services, including: DICOM Receive Storage Commit, and Modality Worklist
- DICOM send
- View DICOM header information
- DICOM CD/DVD import and archival
- DICOM CD/DVD burn functionality
- Ability to edit patient/study information
- Easy to filter, search and sort modality worklist
- Worklist search
- Lossless image compression
- Touchscreen friendly PAC
- Multiple hardware manufacture integration Dual-monitor friendly
- Customizable Report generation User customizable filters
- View images from any web-enabled PC with TigerView Cloud Multiple image viewing

### TigerView Cloud

- Remote viewing
- Disaster recovery
- Multi-office and multi-provider sync
- Off-site Cloud storage/backup
- Centralized database
- Automatic software updates

### PACS Viewer - Modality Support

- Direct Digital Radiography (DR)
- Intra-oral cameras
- Computed Radiography (CR)
- Extra-oral cameras
- Digitized x-ray film
- Color digital camera images for morphing and wound care progression
- Ultrasound (US)

### Measurement & Annotation Tools

- Line Measurement
- Area Measurement
- Angle Measurement
- L/R Cervical Normal Arc
- L/R Thoracic Normal Arc
- L/R Lumbar Normal Arc
- Jackson's Angle
- Baseline Measurement
- Perpendicular Baseline Measurement
- George's Line
- Cobb's Angle
- Stitching
- Automatic Calibration Tool
- Cropping
- Magnification and Zooming
- Window Leveling
- Colorized Density
- Image Negation
- Image Sharpening
- Image Marking and Text Box Features

### EMR or Charting Bridging - Easy link from EMR

- Proprietary bridges
- HL7
- DICOM modality worklists
- TigerLink



For complete product information, visit [televeresystems.com](http://televeresystems.com)