



VIEWWORKS

Choice of digital radiography with Veterinary specific workflow and features

- Wired or wireless flat panel detectors
- 17 x 17 in (43 x 43 cm) / 14 x 17 in (35 x 43 cm) / 10 x 12 in (24 x 30 cm)
- Robust detectors
- Excellent image quality
- Optimized veterinary workflow maximum efficiency
- Fully processed images in ~5 seconds

DEL MEDICAL



#9000-VETFLYER

Rev. A

Digital Radiography

STANDARD UNIVERSAL VETERINARY RADIOGRAPHIC SYSTEM

TABLE

- Stationary Table Top
- Conveniently placed exposure switch and table top release strips are foot activated to ensure safety for the veterinarian and animals being examined

TUBE STAND

- Articulating arm on side of column for mounting generator control console

GENERATOR (single or three phase)

- 30 kW, high frequency, single phase, 400 mA, 125 kVp, with anatomical programming console
- 32 kW, high frequency, three phase, 500 mA, 125 kVp, with anatomical programming console

X-RAY TUBE

- 140 kHU, 125 kVp, 1.0-2.0 mm focal spots

SYSTEM OPTIONS

TABLE

- 2-way float top
- 4-way float top

TUBE STAND

- Fixed SID, 40 in (101 cm) to image receptor
- Variable SID, 41 - 9 in (104 - 23 cm) vertical range with electric locks
- Additional articulating arm to mount DR workstation monitor

GENERATOR (single or three phase)

- 40 kW and 50 kW generator options available upon request
- Pedestal stand or wall mount for control console

X-RAY TUBE

- 200 kHU, 125 kVp, 0.6-1.5 mm focal spots

New York

28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Illinois

241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

www.delmedical.com