

Medlink CDR

Compact, Advanced X-Ray System Designed for
Chiropractic Care That Includes New 45° Arm Rotation



MEDLINK IMAGING

a **VIEWWORKS** company


Compact Chiropractic X-Ray System with 45° Arm Rotation

Medlink CDR, Advanced X-Ray Technology

Medlink has a new upgraded compact x-ray system designed specifically for chiropractic care. The Medlink CDR has an ultra small footprint and comes with a variety of features such as full software integration, **45° arm rotation**, magnetic lock brakes, 40kW generator, dual LED laser collimator and a removable grid.

This x-ray solution can fit in virtually any x-ray room as well as drastically cut down installation time reducing clinical downtime.

1-Year Warranty
(extended warranty available)



The VXVUE acquisition software is an award-winning intelligent software solution which utilizes the most advanced imaging processing. VXvue's image viewer program contains specialized modules to acquire, adjust, and manage diagnostic images created by the VIVIX-S Series.



Bucky Enclosure
Features fully counterbalanced powerful magnetic lock brakes from the handle control to raise and lower the system effortlessly. The bucky also rotates to maintain tube alignment during arm angulation positioning.

MEDLINK IMAGING

a viewworks company

Raytech Diagnostics

<https://raytechdiagnostics.com/>

Tel: 888-721-1982

sales@raytechdiagnostics.com