

# Porta 100HF for Veterinary

“The Green Machine - Dual 100HF”  
High Frequency Portable X-ray Unit

**TELEVERE**  
SYSTEMS



**100 kVp, 30MA**

**Dual Distance Lasers 26" SID**

**Weighs Only 19 lbs  
(Lightest 100 kVp on the market)**

**Six (6) Year Parts Warranty**

**Perfect X-ray Machine for All  
Veterinarians: Small Animal,  
Equine and Mixed Animal**



# Porta 100HF for Veterinary

## “The Green Machine”

### High Frequency Portable X-ray Unit



#### Specifications

“The Green Machine™” Porta 100HF is the newest technology to hit the Veterinary market. It is the only portable x-ray machine on the market boasting a (6) six year parts warranty. Manufactured in Japan with the highest quality parts available. This powerful and light x-ray machine weighs only 19 lbs.

“The Green Machine™” is versatile and perfect for small animal, equine, and mixed animal practices.

#### System Features

- High Frequency – High Output 2 kW, 100 kVp, 30 mA
- Rotary Switches for easy and quick selection of kVp and mAs
- Programmable APR stores up to (8) presets
- Tape measure
- DR Flat Panel compatibility with modular sync connector
- Automatic mAs control for line voltage compensation (helpful with long extension cords)
- Dual Distance lasers preset to 26” providing accurate FFD for Equine Usage
- Aluminum Carrying Case for easy transport
- (6) Year Standard Parts Warranty

#### Product Specifications

- kVp setting: 40-100 kV (in 2 kVp step)  
Rotary switch for digital display with 7 segment LED
- mAs selection: 0.3-20 mAs (32 settings)
- X-ray tube:
  - a) Focal spot: 1.2 x 1.2 mm
  - b) Anode heat storage capacity: 20 kHU
  - c) Target angle: 16 degrees
- Power requirement: 100-120 V or 200-240 V
- Weight: 19.4 lbs

#### Upgrade Options

- Maximal View X-ray Table
- Digital DR Systems – Equine
- Digital DR Systems – Small Animal
- Digital CR Systems

RAYTECH DIAGNOSTICS  
<https://www.raytechdiagnostics.com/>  
[sales@raytechdiagnostics.com](mailto:sales@raytechdiagnostics.com)  
613-799-4171