



ImageWorksVeterinary

Image-Vet FleXX™

FleXXibility to take dental radiographs anytime...anywhere!



Portable Intraoral Dental X-Ray System

ImageWorks Veterinary presents the latest generation in the series of Image-Vet™ X-Rays. Featuring complete portability, the Image-Vet FleXX™ sets the standard for veterinary dental radiography.

Quick and Efficient

Your time is valuable. The Image-Vet FleXX™ allows a complete series of dental radiographs in 65% less time* compared to conventional generators.

Safe and Powerful

At 60kv/2mA, the Image-Vet FleXX™ has the power to produce high quality radiographs on even the largest patients. And it's FDA approved for handheld use.

PRODUCT HIGHLIGHTS

- + Portable - FleXXible to use anywhere
- + Compact and light - only 4 lbs
- + Safe and easy operation with one hand
- + Take more radiographs, faster
- + No assembly, no installation
- + FDA approved for hand-held use
- + Allows radiograph production with one hand
- + No arm to drift
- + Oil-cooled Toshiba tube



ImageWorksVeterinary

Image-Vet FleXX™

FleXXibility to take dental radiographs anytime...anywhere!

"A dental digital radiography machine is the most profitable piece of equipment in a clinic. There are few patients that don't need dental radiographs. And I can't imagine doing dentistry without them."


- DR. KATE KNUDSON (VETERINARY ECONOMICS - Case Studies: Make sure new equipment earns its keep)



The neck strap protects the Image-Vet FleXX™ and allows for simple and easy one-handed operation, allowing the user total FleXXibility to operate acquisition software, position the digital sensor or tend to the patient.



The Image-Vet FleXX™ requires NO installation or assembly. Its two rechargeable batteries, neck strap and back scatter shield all fit conveniently in the supplied hard-shelled transport case (not pictured).



The Image-Vet FleXX™ integrates seamlessly with EVA Vet Digital Dental Imaging Systems and can be used with other higher quality digital sensors and high speed dental films. Ask about a package deal for greater savings!


PRODUCT SPECIFICATIONS

- High Frequency 60kV / 2mA
- Compact & Light - only 4 lbs
- Double lead shielding prevents x-ray leakage
- High capacity Lithium-ion battery (up to 500 shots per charge)
- Oil-cooled Toshiba tube
- FDA approved for hand-held use
- Generates less harmful radiation than conventional systems

ACCESSORIES

- Neck strap for safety and FleXXibility
- Back scatter shield for user safety
- Two rechargeable battery packs
- Battery charging station
- Hard-shelled transport case included

AUTHORIZED DISTRIBUTOR



RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
 613-799-4171